







Instructions for Use




Formfit® Tracker

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	7
FR	Notice d'utilisation	9
ES	Instrucciones para el uso	12
IT	Istruzioni per l'uso	14
NO	Bruksanvisning	16
DA	Brugsanvisning	18
SV	Bruksanvisning	20
FI	Käyttöohjeet	23
NL	Gebruiksaanwijzing	25
PL	Instrukcja użytkowania	27
TR	Kullanım Talimatları	29
JA	取扱説明書	32
ZH	中文说明书	34
KO	사용 설명서	36

SYMBOLS

Medical Device

LIST OF PARTS

Fig. 1:

- a. PowerLock™ strap
- b. Hook and loop surfaces
- c. Strap ends
- d. ReflexWing™ buttress
- e. Kneecap opening
- f. CustomFit™ hinges
- g. Hinge pocket

INTENDED USE

The device is intended for external support or stabilization of the patella

The device must be fitted and adjusted by a healthcare professional.

Indications for use

For knee conditions that may benefit from realignment of the patella and/or patients that would benefit from an added lateral support of the patella, such as:

- lateral patellar subluxation or dislocation
- lateral maltracking

No known contraindications.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

FITTING INSTRUCTIONS

Fastening dots

Dots show you where to fasten the strap ends. Fasten parts of the device with two vertical dots in the same area. Fasten parts with single dots in the same area.

Note: Make sure to fasten the strap ends when the device is not in use. If you do not fasten the strap ends, they can damage the device.

Device Application

1. Hold the device by the hinge pockets and pull it up over your knee (**Fig. 2**).
2. Make sure your kneecap is in the middle of the opening (**Fig. 3**).
3. Use fingers to keep the buttress pushed against the side of your kneecap (**Fig. 4**). Pull the bottom strap around the back of your knee to tighten it first. Then fasten the bottom strap as the dots show you.

Caution: Do not tighten too much. If you tighten too much, you can cut off your blood flow.

4. Repeat for the top strap (**Fig. 5**).
5. The buttress is now in position and pushes tightly against the side of your kneecap.

Device Removal

1. Pull the strap ends off of the device and fasten them to the strap (**Fig. 6**). Fasten as the dots show you.
2. Hold the brace by the hinge pockets. Then pull it down to take it off.

Device Adjustments

Using the buttress

If necessary, you can bend and move the buttress to adjust it to your patella (**Fig. 7**).

Moving the buttress

1. Pull the buttress off of the loop material (**Fig. 7a**).
2. Put it in the correct position for your patella.

Bending the buttress

1. Pull the strap ends off of the device. Then fasten the strap ends as the dots show you.
2. Fully open the device.
3. Push the middle of the buttress down with your index finger (**Fig. 7b**).
4. Lift the end you want to adjust.
5. Bend it to agree with the shape of your patella.
6. Push the end back down again.

Shortening the straps

If the straps are too long, you can make them shorter.

Note: Only shorten the strap in single steps of 0,5 cm/¼". If you cut off too much, you may not be able to use the brace correctly.

1. Fully remove the strap end from the strap.
2. Then use scissors to cut off a short piece of the strap (**Fig. 8**).
3. Put the device on again to check how it fits. If the strap is too long, repeat steps 1) and 2). If the strap length is good, move to step 4).
4. Put the strap end back on the strap again.

Using the hinges

Removing the hinges

If necessary, the hinges can be heated and bent by a certified prosthetist orthotist (CPO) to agree with the shape of your leg. To remove the hinges, do as follows:

1. Pull the strap ends off of the device and fasten them to the strap (**Fig. 6**). Fasten as the dots show you.
2. Push the hinges out of the device through the openings on the hinge pockets (**Fig. 9**).

Putting the hinges back

To put the hinges back into the device, do as follows:

1. Pull the strap ends off of the device and fasten them to the strap (**Fig. 6**). Fasten as the dots show you.
2. Push the hinges into the device through the openings on the hinge pockets (**Fig. 9**).

USAGE

Cleaning and care

- Machine-wash with a mild detergent at cold temperature (30°C).
- Air dry.

Note: Do not tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

Note: Secure the hook and loop fasteners before washing to avoid damaging the textiles.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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