


ENGLISH

MD Medical Device

INTENDED USE

The device is intended for unicompartmental unloading of the knee. The device must be prescribed by a healthcare professional.

Indications for use

- Unicompartmental osteoarthritis of the knee.
- Other conditions that may benefit from unicompartmental unloading, such as degenerative meniscal tears.

No known contraindications.

Warnings and Cautions:

- Regular healthcare professional supervision is recommended for patients with peripheral vascular disease, neuropathy, and sensitive skin.
- Make sure the device fits properly to minimize the possibility of skin irritation. Increase usage time gradually as the skin adapts to the device. If redness appears, temporarily decrease the usage time until it has subsided.
- Care should be taken not to overtighten the device.
- Make sure the device fits properly to achieve effective pain relief.
- The use of the device may increase the risk of deep vein thrombosis and pulmonary embolism.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device. Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, excessive pressure, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

FITTING INSTRUCTIONS

While carrying out the following instructions, please refer to the overview figure for locating components mentioned in the text (**Fig. 1**).

Device Application

- Ask the patient to sit down on a chair and partially bend the knee with the foot flat on the floor. Position the fully opened device on the leg so that the Hinge (A) is on the affected side of the knee and the Kneecap Opening (B) is positioned around the kneecap (**Fig. 2**).

- Ensure proper alignment of the device on the leg (**Fig. 3**).
 - Height positioning: Lower edge of the Kneecap Opening should be below the kneecap (3a) and align Össur logo on the Hinge with upper edge of the kneecap (3b).
 - Side positioning (3c): Front edge of the Hinge cap should be at midline of leg.
- Attaching the Tabs (C, D, E, F).
 - While keeping the device in place with

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from

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